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rev03 04-12-2025

BIREADER

Incubator for metrics Rapid Biologicals

Product Overview

The BI READER TecSpores® incubator is engineered for the automated reading of TecSpores® Rapid Biological Indicators utilized in the validation of steam sterilization (STEAM), ethylene oxide (EO), vaporized hydrogen peroxide (VH₂O₂), and low-temperature formaldehyde (LTSF) processes.

Characteristics

- Automatic fluorescence analysis of Rapid TecSpores® IB.
- Intuitive touch interface yielding results in five steps.
- Swift and dependable identification of process efficiency.
- Compatible with various sterilization protocols: EO, STEAM, VH₂O₂, and LTSF.
- Visible results on screen with clear indications: (+) Inadequate sterilization / (-) Effective sterilization.
- It necessitates the utilization of a rapid, unprocessed IB as a positive control.

Specs:

- Incubation temperature:
- 37°C (EO - 240 minutes)
- 58 °C (VH₂O₂ / STEAM - 20 minutes, LTSF - 60 minutes, STEAM - 180 minutes)
- Compatibility: IB Rapids TecSpores®

How to Utilize:

Activate one processed and one unprocessed IB using the breaker. Position the IB in the incubator, select the channel and program. Await the result on the screen: “+” signifies failure, while “-” denotes success. The unprocessed IB should exhibit “+” as a control.

For the sole use of healthcare professionals and institutions.

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