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Sterilization indicator on Trio 2080 adhesive label

Product Overview

EFELAB self-adhesive labels are specifically engineered for use with printing and coding systems, facilitating the identification of materials subjected to ethylene oxide, moist heat, and dry heat sterilization. Their sensitive ink responds to the critical parameters defined in ISO Standard 11140-1.

Features

- Dimensions: 80 mm in length x 40 mm in height.
- Compatible with conventional market coding systems.
- Composed of 60 to 70 g/m² construction paper.
- Reinforced adhesive featuring high initial tack, exceptional adhesion, and medium cohesion.
- Printing with two 4 mm sensitive ink bands on the right (DRY - STEAM) and left (EO) margins.
- Dual-sensitive ink that alters color in response to moist heat (steam) and dry heat (oven) sterilization processes.
- Specialized sensitive ink for ethylene oxide sterilization processes.

Presentation:

They are offered in individual packages (reel or roll) containing 1,000 self-adhesive labels, affixed to a cardboard cone with a diameter of 25 mm and a maximum overall diameter of 95 mm.

Sensitive dye

The sensitive ink transitions to brown when subjected to dry heat at $160^{\circ}\text{C} \pm 5^{\circ}\text{C}$ for a minimum duration of 10 minutes ± 45 seconds. Additionally, it responds to moist heat sterilization processes in steam sterilizers in accordance with ISO 11140-1, operating at a pressure of 1.1 bar at 121°C for 10 minutes ± 5 seconds, or 2.2 bar at 134°C for 2 minutes ± 5 seconds.

The sensitive ink undergoes a color change when subjected to the ethylene oxide sterilization process, responding to the critical parameters defined in ISO 11140-1. These parameters encompass a concentration of $600 \text{ mg/L} \pm 30 \text{ mg/L}$, a temperature of $30^{\circ}\text{C} \pm 1^{\circ}\text{C}$, a minimum exposure duration of 30 minutes ± 15 seconds, and a relative humidity of $60\% \pm 10\%$.

How to Utilize

Ensure that the packages are clean, securely packed, and permit the circulation of the sterilizing agent.

The print on the right margin transitions from greenish-blue to brown following sterilization with either moist or dry heat. The print on the left margin shifts from light blue to yellow/ochre when exposed to ethylene oxide.

If there is no alteration in color, the sterilization process has been unsuccessful.

Precautions

This product should not be utilized for any processes other than moist heat in an autoclave (vapor), dry heat (oven), or Ethylene Oxide.

Storage

Store in a cool environment. Maintain dryness and protect from light.

Useful Life

The product has a shelf life of 60 months from the date of manufacture.

Regulation

All activities conducted at EFELAB adhere to the ISO 13485 quality standard and comply with Good Manufacturing Practices (GMP) to guarantee product quality.

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