INSTRUCTIONS FOR USE EFELAB STERILIZING GAS AMPOULES

Nro. Rev. 02

1-Do not open the plastic wrapping bag in which each individual Efelab ampoule is sealed

Each ampoule has a cover shield around the ampoule. The ampoule and shield are sealed in a plastic bag. The plastic bag is a diffusion membrane of known permeability whose function is to contain the gas given off by the ampoule and release it at a controlled rate during the sterilization cycle. The cover shield around the ampoule prevents broken glass from the opened ampoule from puncturing the gas releasing bag. Efelab sterilizing gas allows to sterilize thermo labile materials, accounting for products that cannot be submitted to either steam or dry heat sterilization processes. Efelab sterilizing gas penetrates into and easily spreads throughout different materials, such as plastics, rubber, artificial fibers, textiles, papers, among others.

2-Preparation of material for sterilization

Material to be sterilized by EFELAB sterilizing gas must be surgically clean and dry. Coating of dry proteins like dry blood, pus or feces protect microorganism and delay sterilization process. Washing contaminated instruments with soap, detergent or another cleaning product, and water prior to exposure to sterilizing gas rehydrates desiccated spores and is strongly recommended as a routine procedure. If the nature of materials is such that the aforementioned water treatments are harmful, then a pretreatment of at least one hour in a chamber having saturated humidity will be necessary for an overall pretreatment routine.

Although the active ingredient contained in the ampoule is a highly diffusible gaseous sterilizing agent, occlusive caps, plugs or stylets should be removed from instruments so that the gas can freely penetrate. Do not pack the container so tightly with cloth or gauze that the gas diffusion is compromised. Hollow-bore needles and plastic or rubber tubing must be open and free from plugs. Syringes should be packaged disassembled with their plungers out of the barrel. Water on instruments at the time of exposure to sterilizing gas may react with it and reduce its effectiveness. Please note that the material to be sterilized must be perfectly clean and dry before wrapping and processing. Place each element to be sterilized into special sterilization packaging to avoid their contamination after the sterilization process. Efelab polyethylene film with built-in sterilization indicator and medical reels or pouches are recommended as packaging material. Close each package with a proper wrapping material to avoid any contamination after the sterilization process.

3-Sterilization method

Have the material to be sterilized prepared as described above and ready to place in the liner bag. Take out a liner bag from the container box and place it fully opened within a sterilizer or a special container. A new liner bag must be used with each load. Get one sterilizing gas ampoule and unroll the plastic bag. Do not open or break the plastic bag in which it is wrapped, - like it was mentioned in step 1 (one) - because this is a gas-diffusing membrane. The ampoule is prescored in its neck to facilitate its opening and the following gas release. Grab the ampoule, identify its shape through the bag and cover shield and snap off the neck of the ampoule. This releases the gas within the ampoule's bag. Place the ampoule, still wrapped in its sealed plastic bag,

inside the bottom of the liner bag. Place the material to be sterilized on top of the gas-release ampoule bag. Once the bag was filled with material twist the liner bag shut and hold tight the seal with one plastic seal which is also provided with this sterilization kit.

Place the container bag into a special container, sterilizing cabinet or chamber and record the date and hour the sterilization has begun. Put the container in a safe and well-ventilated place, away from flammable items or sparks. Leave the container at room temperature -at least 20°C or 68°F- for 12 hours before removing the sterilized material. Once sterilization cycle is completed, open the container bag and remove the sterile material. Both container bag and other elements used (the exhausted and empty glass ampoule and its bag) may be disposed off in an ordinary rubbish dump.

4-Special materials

When sterilizing long length plastic or rubber tubing - e.g. urethral catheters, heart catheters, IV extension tubes - use two glass ampoules and increase the sterilization cycle time to 24 hours.

5-Special considerations

Conditions to take into consideration when sterilizing with EFELAB ethylene oxide gas ampoules without a sterilizer machine:
Room temperature has to be at least 20° C or 68° F.
Relative humidity should not be less than 30 %.
Contact time with sterilizing agent -sterilization cycle- shall last 12 hours or more.

Personnel should wear personal protection equipment during these procedures.

6-Precautions

Danger! Keep out of reach of children. For professional use only.

Efelab sterilizing gas is extremely flammable, keep away from heat, sparks and flame. Do not smoke near the sterilization container or machine while loading or unloading it. Remove batteries from electrical instruments being sterilized and wrap them separately to avoid the possibility of an electrical spark igniting the gas during the sterilization cycle. Keep sterilization container closed. Avoid breathing Efelab sterilization gas vapors in concentrated amounts. Efelab sterilizing gas is as irritating to the lungs and mucus membrane as is ammoniac gas. As with other chemical vapors there is a change of an occasional allergic response to this product in a sensitive individual. Such individual should not handle this product nor breather its vapors nor allow materials sterilized in it to come in contact with their skin, eyes and clothing. All users should avoid contact with skin, eyes and clothing, including shoes. Flush skin or eyes with plenty of water for at least fifteen minutes in case of accidental contact and get medical attention. Failure to adequately aerate gas absorbing materials may lead to contact chemical burns. Foam rubber, plastic foams ad sealed air cushioned devices like some anesthesia mask must be aerated at least 24 hours before use. Sterilizing agent penetrates and adheres to rubber and PVC and it is progressively released after the sterilization cycle. Hence all the material sterilized with this gas shall enter into quarantine time period and it should be efficiently

ventilated. According to the nature of the product, ventilation of sterilized material shall last 24 hours or more. Take into consideration that sterilized materials composed of porous or high-absorbent elements must be aerated for at least 24 hours. If sterilized materials are going to be used within or over the human body, they have to be ventilated for at least 72 hours. Ventilation time may be substantially reduced by using a special heated aerator or a special ventilation cabinet or a proper sterilization machine. For aeration times, using a heated aerator or a sterilizer machine, please consult the aerator or machine users' manual. It is not necessary to remove such material from cloth or paper wrappings to obtain adequate diffusion of the residual. Vacuum dehydration, chemical desiccation or prolonged exposure to low relative humidity -below 30 %- has been demonstrated to produce spores highly resistant to sterilization. Rehydration to spores so changed and hence reversion to normal sensitivity does not seem to occur until they have been actually wetted or placed in essentially 100% relative humidity. Do not attempt to sterilize materials which may be carrying dry spores without first soaking them in warm water or washing the material with water and soap or detergent. If the nature of material is such that water treatment specified above is harmful, then a pretreatment in a chamber having saturated humidity will be necessary for an overall pretreatment routine of at least 60 minutes.

Efelab sterilizing gas is a powerful polymerizing agent; it should not be used for sterilizing food-stuffs or drugs. Sterilization efficacy is extremely important. Hence it is strongly recommended that the user establish a routine for monitoring each sterilization cycle. Should you require further information, please refer to the MSDS of the product or visit our FAQ's section at www.efelab.com.

7-Important: Variations from these Instructions

Any variation from these instructions is done at the user's risk. Efelab sterilizing ampoule is designed to be used only in the liner bag provided with the ampoule and only in proper sterilization container or sterilizer machine.

8-Technical description of the system

Appropriate sterilization containers and machines serve as a guard against inadvertent ignition of the contents by sparks or open flames during the sterilization cycle. The liner bag, when properly twist closed has a capacity of 11 liters, 20 liters or 34 liters depending on the ampoule gas content: 5 cc, 10 cc or 20 cc respectively. It serves as a gas diffusion membrane, retaining sterilizing gas enough time to sterilize its content and then diffusing it into the surrounding container or room air at a slow rate. Each ampoule releases approximately 4.495 mg (Product Code 4005), 8.890 mg (Product Code 4010) and 17.890 mg (Product Code 4020) of sterilizing gas at room temperature (20°C or 68°F) at sea level atmospheric pressure giving a minimum peak concentration within the liner bag of 500 mg/1000cc. Laboratory tests confirm that this configuration will kill most resistant spores known (provided they have not been desiccated) with a 12 hours contact-time-sterilization-cycle..

9-Storage, Shelf life

Efelab sterilizing gas ampoules must be stored in a cool place and protected from direct sunlight. Under normal conditions its shelf life will exceed two

years. As long as the material contained into the ampoule is liquid at 20°C or 68°F, it is powerful enough to use. Store this product in accordance with national or local storage procedures. For further information contact your national or local authorities.

10-Product presentations

Product Code	Ampoule Content		Container Bag		Presentation
	Volume	% ETO	Measurement	Capacity	
4005	5 cm ³	90%	31 x 60 cm	11 liters	20 units per box
4010	10 cm ³	94 %	45 x 78 cm	20 Liters	10 units per box
4020	20 cm ³	97%	56 x 96 cm	34 liters	10 units per box

E.M. DE SALUD Y A.S. 21129-93-5
MADE IN ARGENTINA

For further information contact us at:

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