



Indicator Tapes for Hydrogen Peroxide sterilization (VH2O2).

Cellulose free substratum tape with hi-performance chemical indicator ink printed that shows a noticeable change of color after being exposed to Hydrogen Peroxide plasma sterilization conditions, thus indicating the operator of CSSD's that the material has been effectively processed.



Main features of this product:

- ◆ Hi-quality indicator ink that shows a very noticeable color change after being processed.
- ◆ Easy to use in order to mark and fasten sterilization wrapping materials.
- ◆ Manufactured with synthetic compounds (cellulose free).
- ◆ Adhesive compatible different types of wrapping materials.
- ◆ Every tape has a label in its inner cone showing the sterilization type the product was developed for and the color the ink shall turn to after being submitted to a Hydrogen peroxide plasma sterilization cycle, which avoids confusion.