



## Sterilization Indicator in Self-Adhesive Dot Label

Self-adhesive labels with chemical sterilization indicator are used to identify different packages or materials that will be exposed to a sterilization process. Show a color change noticeable and defined after being processed under suitable conditions sterilization. This circular model also supports the possibility of using a labeling machine to apply each label makes it a quick and easy method to identify materials..


Code / Catalog No.	Size	Sterilization Type:
2001	15mm / 0.59in	Ethylene Oxide



Before



After



### Main Features:

- ◆ Adhesive paper that provides good adhesiveness, even at high temperatures.
- ◆ High fidelity indicator Inks. Complies with ISO 11140 / Type I.
- ◆ Labels printed with information related to the sterilization process for which the product it has been developed for and the target color the ink will show after being exposed.
- ◆ Container boxes printed with product lot and expire date for a correct product traceability.
- ◆ 1000 labels per roll.