



## Continuous rotary sealer for medical reels



Continuous sealing machine of rotary type developed for closing medical reels, heat sealable wrappings of medical grade and pouches of paper of Tyvek®

### **Main features of this device:**

Body made of rigid anodized aluminum and top cover of AISI 304 stainless steel, materials that minimize its weight but provide maximum robustness, it is one of the most resistant and compact models in the market

Automatic start system by material detection at the beginning that maximizes the durability and conservation of all the mechanical components of the equipment

Precision electronic controller that allows optimal regulation and adjustment of the sealing temperature, with digital display.

Internal parts manufactured in CNC and bearings of first quality.

Approximate working speed of 9 meters per minute.

Lateral slide with rule for a precise adjustment of the margin to edge of the seal.

Sealing band of 3 lines of 12 mm width.

Total power: 550 Watts, AC 220V, 50 Hz. Working range from 50°C to 220°C.

Measurements: 435 x 235 x 135 mm. Weight 8.6 kg.