




EFELAB



Quality in Sterilization Supplies



ETHYLENE OXIDE CARTRIDGES

EO GAS IS A WELL-KNOWN AND POWERFUL BACTERICIDAL, VIRICIDAL Y FUNGICIDAL AGENT.



Main Features:

- ◆ High purity cartridges, contain 100% ethylene oxide.
- ◆ They are hermetic and inviolable, ISO 13485 certified.
- ◆ They do not contain CFC (Freon), hence they do not affect the ozone layer.
- ◆ They are safe for operators working at the plant.
- ◆ Manufactured under strict quality controls, verifying standard, weight and technical characteristics of the product, on unit by unit bases, being one of the safest alternatives available in the market.
- ◆ Available in several sizes, there is a suitable choice for almost every sterilizer available in the market.

CODE	MEASUREMENTS		CONTENT
PRODUCT	DIAMETER	HEIGHT	GRAMS
1007/10	25 mm.	85 mm.	7/10 grs.
1030	25 mm.	105 mm.	30 grs.
1067	35 mm.	105 mm.	67 grs.
1070	35 mm.	140 mm.	70 grs.
1100/3	38 mm.	165 mm.	100 grs.
1128/3	38 mm.	165 mm.	128 grs.
1130	50 mm.	133 mm.	130 grs.
1134/3	38 mm.	165 mm.	134 grs.
1170-3	45 mm.	176 mm.	170 grs.
1180	50 mm.	174 mm.	180 grs.
1200	50 mm.	174 mm.	200 grs.
1230	50 mm.	174 mm.	230 grs.
1250	65 mm.	260 mm.	250 grs.

