



See-Through Flat Reels



//Blue Peel//
STERILIZATION PACKAGING

Manufactured with 60 GR/M2 medical grade paper, which has been specifically developed for the production of peeling-reels under European norms, and bi-laminated transparent film with anti-bacteria barrier. This conjunct provides optimum resistance, durability y peeling.

- ◆ Built-in chemicals indicator for Steam and Ethylene Oxide (EO) processes.
- ◆ Easy visualization of the packaged contents though its bi-laminated film.
- ◆ Presentation: Gusseted and flat reels.
- ◆ Manufactured under UNE-EN 868 Norm; Part 5 of the European Community (CE).

Main characteristics of the medical grade paper used to manufacture Blue Peel reels.

- ◆ Manufactured with pure paper cellulose, chloride and optic whiteners free, it is environmentally friendly.
- ◆ Fulfill the requirements of the strict CE norm; EN 868 (parts 1,2,3,4 and 6) – specific for medical grade papers.
- ◆ High performance anti-bacteria barrier: it provides more security against contamination of the package.
- ◆ High mechanical resistance: It provides security when it comes to materials with acute edges or heavy objects.
- ◆ Pore size controlled: Necessary condition to assure the maintenance of sterility inside the package.
- ◆ High resistance against liquids which provides more security against contamination
- ◆ Superior sealing/peeling quality.

