



Polyethylene Film

with built-in sterilization Indicator for Ethylene Oxide (EO).



EFELAB polyethylene film accelerates the preparation of materials to be processed as it has a highly reliable and durable built-in sterilization indicator that complies with ISO 11140 Standards. Its colour change is very significant -*it turns from violet to light green*. Because of the polyethylene transparent-appearance, the materials to be sterilized may be easily displayed and handled. Once processed, this film maintains its color even in presence of sunlight, thus allowing correct material identification. It can be used as a sterilization record.

Main Features:

- ◆ As it is a transparent packaging it allows to easily see the contained materials.
- ◆ It provides sterilized materials with a long-lasting maintenance period.
- ◆ It allows the operator to properly adjust the packaging length.
- ◆ It is an inexpensive alternative compared to other sterilization wrapping materials.
- ◆ It comes in rolls of 60 mm; 80 mm, 100mm; 15 mm; 20 mm and. 300 mm width; 100 meters long and 60 microns thick.