



Sterilization Indicator Tapes

Indicator tapes are used to fasten packaging, container and different sterilization wrapping materials. They can also be used to close envelopes, see-through flat reels and packages when a sealing machine it is not available. Printed with a high quality chemical indicator ink they show a notorious colour change after being exposed to proper sterilization conditions. Therefore show the operator that the materials have been processed. Available for different sterilization processes: EO gas, Steam and Dry Heat.



Code/ Catalog No.	Measurements	Sterilization Process:
3012	50 mts. X 12 mm.	Steam
3012 Y	60 yardas. X 1/2 in	Steam
3015	50 mts. X 15 mm.	Steam
3019	50 mts. X 15 mm.	Ethylene Oxide (EO)
3018	50 mts. X 15 mm.	Dry Heat (DRY)
3017	50 mts. X 18 mm.	Steam
3017 Y	60 yardas. X 3/4 in	Steam
3016	50 mts. X 18 mm.	Ethylene Oxide (EO)
3020	50 mts. X 18 mm.	Dry Heat (DRY)
3025	50 mts. X 24 mm.	Steam
3025 Y	60 yardas. X 1 in	Steam

Main Features:

- ◆ Latex-treated Crepe Papel.
- ◆ Rubber-based special adhesive that provides high adhesiveness at high temperatures.
- ◆ Lacquer on surface to facilitate unrolling.
- ◆ High quality indicator inks.
- ◆ Every tape is printed in the inner cone with type of sterilization process to be used with and the colour the indicator ink turns to, thus avoiding misunderstandings.
- ◆ Container bags printed with lot number and validity for correct product traceability.